



SOLUTIONS

PRODUCTS

RESOURCES

EVENTS

CONTACT US

REQUEST QUOTE

# Understanding Dynamic Light Scattering

When in solution, macromolecules are buffeted by the solvent molecules. This leads to a random motion of the molecules called Brownian motion. For example, consider this movie of 2  $\mu\text{m}$  diameter particles in pure water. As can be seen, each particle is constantly moving, and its motion is uncorrelated with the other particles. (Movie courtesy of Dr. Eric R. Weeks, Emory University).



Brownian motion.

As light scatters from the moving macromolecules, this motion imparts a randomness to the phase of the scattered light, such that when the scattered light from two or more particles is added together, there will be a changing constructive or destructive interference. This leads to time-dependent fluctuations in the intensity of the scattered light (Fig. 1).

In dynamic light scattering (DLS), the time-dependent fluctuations in the scattered light are measured by a single photon counting module. The rate of fluctuations is directly related to the rate of diffusion of the molecule through the solvent, which is related in turn to the particles' hydrodynamic radii. Smaller particles diffuse faster, causing more rapid fluctuations in the intensity than larger particles. Therefore, the fluctuation in light intensity contains information about the diffusion of the molecules and can be used to extract a diffusion coefficient and calculate a particle size.

DLS is employed by the DynaPro® NanoStar®, the DynaPro® Plate Reader, the Mobius™ and the WyattQELS™ module for MALS detectors to determine the effective particle size. The analyte's translational diffusion coefficient,  $D_t$ , is obtained by automated nonlinear least squares fitting of the autocorrelation function that quantitatively describes the measured time-dependent fluctuations in light scattering intensity. The analysis is done directly in the accompanying DYNAMICS®, DYNAMICS® Touch™ or ASTRA® software. A typical autocorrelation function for a monodisperse sample is shown in Fig. 2.

The Stokes-Einstein equation then gives the hydrodynamic radius,  $R_h$ , (Fig. 3) corresponding to the measured  $D_t$ :




Figure 1. Brownian motion results in measurable fluctuations in the light intensity over time.




Figure 2. Autocorrelation function for a monodisperse sample with 50 nm  $R_h$  particles.

$$R_h = \frac{kT}{6\pi\eta D_t}$$

where  $k$  is Boltzmann's constant,  $T$  is the temperature in K, and  $\eta$  is the solvent viscosity.

Different sizing techniques, e.g., DLS, small angle X-ray scattering, microscopy, and molecular modeling may report different types of radii. It is therefore important to know how a reported "size" was determined and whether it refers to the radius or diameter of the molecule. The  $R_h$  measured by DLS is the radius of a hard sphere with the same diffusion coefficient as the sample. Other measures include  $R_g$  (the radius of gyration, or root-mean-square radius, obtained by e.g., MALS),  $R_{Vol}$  (the radius of a hypothetical sphere that occupies the same volume as the macromolecule), and  $R_{Rot}$  (the radius subtended by rotating the macromolecule). Fig. 3 illustrates how these different measures compare for the compact, globular protein lysozyme.




Figure 3. (Left)  $R_h$  is the radius of a hard sphere with the same diffusion coefficient as the sample. (Right) Compact molecules, e.g., globular proteins such as lysozyme, generally have  $R_g < R_h$ .

- + Analyzing Particle Size Distributions
- + Measuring Nanoparticle Concentration
- + Characterizing Intermolecular Interactions
- Understanding the Autocorrelation Function

In DLS, the fluctuations in light intensity measured over time are quantified via a second order correlation function  $g^{(2)}(\tau)$ . The function of intensity is shifted by a delay time ( $\tau$ ) and the autocorrelation function  $g(\tau)$  is calculated. Qualitatively, the autocorrelation function is a measure of how similar the intensity function is to itself when shifted by time  $\tau$ . As the value of  $\tau$  increases, the function reaches a baseline of 1.

More specifically, and as described in various light scattering texts (cf. Chu, B. Laser Light Scattering: Basic Principles and Practice; Academic Press: Boston, 1991), the correlation function for a monodisperse sample can be analyzed by the equation:

$$g^{(2)}(\tau) = 1 + \beta \exp(-2\Gamma\tau)$$

where  $\beta$  is the correlation function amplitude at zero delay,  $\Gamma$  is the decay rate, and the baseline of the correlation function relaxes to a value of 1 at infinite delay.

DYNAMICS and ASTRA use a nonlinear least squares fitting algorithm to fit the measured correlation function to equation 2 to retrieve the correlation function decay rate  $\Gamma$ . From this point,  $\Gamma$  can be converted to the translational diffusion coefficient  $D_t$  for the particle via the relation:

$$D_t = \frac{\Gamma}{q^2}$$

Here,  $q$  is the magnitude of the scattering vector, and is given by:

$$q = \frac{4\pi n_0}{\lambda} \sin\left(\frac{\theta}{2}\right)$$

where  $n_0$  is the solvent index of refraction,  $\lambda_0$  is the vacuum wavelength of the incident light, and  $q$  is the scattering angle.

Finally, the diffusion coefficient can be interpreted as the hydrodynamic radius  $R_h$  of a diffusing sphere via the Stokes-Einstein equation:

$$R_h = \frac{kT}{6\pi\eta D_t}$$

where  $k$  is Boltzmann's constant,  $T$  is the temperature in K, and  $\eta$  is the solvent viscosity (Fig. 3).

[View DLS Webinars](#)

[View DLS Application Notes](#)

## RESOURCES

- [Store](#)
- [Email Signup](#)
- [Bibliography](#)
- [Contact Us](#)
- [Careers](#)
- [Blog](#)
- [Request Product Info](#)
- [Webinars](#)
- [Support Login](#)
- [Light Scattering University](#)

Wyatt Technology is the recognized leader in light scattering instrumentation and software for determining the absolute molar mass, size, charge and interactions of macromolecules and nanoparticles in solution.

Wyatt's line of multi-angle static light scattering products couple to size exclusion chromatography (SEC-MALS), field-flow fractionation (FFF-MALS) and stop-flow composition-gradient systems (CG-MALS). Our dynamic light scattering (DLS) products operate in traditional cuvette as well as on-line and automated, high-throughput modes. We also offer unique instruments for electrophoretic light scattering (MP-PALS), differential refractometry and differential viscosity.

+1 805-681-9009 Headquarters | 6330 Hollister Ave, Santa Barbara, CA 93117

[Global Offices](#) | [Site Map](#) | [Feedback](#) | [Terms of Use](#) | [Privacy Policy](#)

© Wyatt Technology Corporation. All rights reserved.

[Privacy & Cookies Policy](#)

Also of Interest

- Introduction to Light Scattering Theory
- Dynamic Light Scattering (DLS)
- Size

Ultimate Guide

Thirteen key questions to ask when buying a Multi-Angle Light Scattering (MALS) or Dynamic Light Scattering (DLS) instrument.

Request Ultimate Guide

Solutions Guide

Our Solutions Guide combines Wyatt's product lines, applications and services into a single comprehensive booklet.

Request Solutions Guide

Biologics Guide

Our Biologics Guide highlights tools and techniques for characterizing biologics and complex drugs.

Request Biologics Guide

RESOURCES

- Store
- Email Signup
- Bibliography
- Contact Us
- Careers
- Blog
- Request Product Info
- Webinars
- Support Login
- Light Scattering University

Wyatt Technology is the recognized leader in light scattering instrumentation and software for determining the absolute molar mass, size, charge and interactions of macromolecules and nanoparticles in solution.

Wyatt's line of multi-angle static light scattering products couple to size exclusion chromatography (SEC-MALS), field-flow fractionation (FFF-MALS) and stop-flow composition-gradient systems (CG-MALS). Our dynamic light scattering (DLS) products operate in traditional cuvette as well as on-line and automated, high-throughput modes. We also offer unique instruments for electrophoretic light scattering (MP-PALS), differential refractometry and differential viscosity.

+1 805-681-9009 Headquarters | 6330 Hollister Ave, Santa Barbara, CA 93117  
Global Offices | Site Map | Feedback | Terms of Use | Privacy Policy

© Wyatt Technology Corporation. All rights reserved.

Privacy & Cookies Policy